


AT10/T2000D/T10D Grip Catalog

Other grip options are available upon request

P1238 56 mm Dia. Aluminum Compression Platen (set of platens)



P1239 116 mm Dia. Aluminum Compression Platen (set of platens)



G0228 O-Ring 15mm/25mm Dia. w/120v option Grips

Other diameter spools available upon request

S1684 – Spool, O-Ring 18.3mm Diameter

S1685 – Spool, O-Ring 22.3mm Diameter



G0227 Eccentric Roller Grips 1kN Aluminum (set of grips)



G0230 Grips, Pneumatic Wedge 10kn Serrated Jaws (set of grips)
Grips come standard with F0405




F0405 - 0-6mm Opening Wedge Pyramid (Serrated) Jaws (supplier p/n TH175-10-BP6) (set of jaws)



G0229 - Grips, Pneumatic Vise 1.2kN Serrated Jaws (set of grips)
Grips come standard with F0403



F0403 - 30x30 mm Vise Pyramid (Serrated) Jaws (supplier p/n TH240k-BP) (set of jaws)
F0404 - 30x30 mm Vise Rubber Jaws (supplier p/n TH240k-BG) (set of jaws)



Filament, Fibers, Yarns, Cords, Wires and Strands Pneumatic Grips (set of grips)
Grips come standard with faces listed below each grip. Additional grip faces are available.

G0231 – Pneumatic Rope 200N Grips
F0406 – 200N Diamond Coated Rope Grip Face

G0232 – Pneumatic Rope 2kN Grips
F0407 – 2kN Serrated Rope Pyramid Faces

G0233 – Pneumatic Rope 5kN Grips
F0408 – 5kN Serrated Rope Pyramid Faces



N105C Pneumatic Closing Grip (set of grips)



G0234 – Grip, EN 1372 Peel (1 grip)

10kN Roller Peel Fixture, roller diameter (20mm), 0 to 45 degree



G0235 – Grip, ASTM D3167 Peel (1 grip)

10kN Roller Peel Fixture, roller diameter (25.4mm), roller width (25.4mm), 0-115 degrees

G0236 – Grip, EN 1464 Peel (1 grip)

10kN Roller Peel Fixture, roller diameter (25mm), roller width (30mm), 0-90 degrees



G0258 – Grips, Pneumatic Vise 3.5Kn 50 x 100 serrated Jaws

G0259 - Grips, Pneumatic Vise 3.5Kn 50 x 100 Wave Jaws

G0260 - Grips, Pneumatic Vise 3.5Kn 50 x 100 Diamond Jaws


G0261 - Grips, Pneumatic Vise 3.5Kn 50 x 100 Rubber Jaws



G0262 – Grips, Pneumatic Vise 2kN 25x50mm Wave Jaw (Wave Jaws included)
Other Grip Faces “Jaws” available upon request



G0263 - Grips, Pneumatic Vise 2.4kN 30x50mm Serrated Jaw (Serrated Jaws included)
Other Grip Faces “Jaws” available upon request



F0457 – Faces, Wave 30x50mm



G0267 – Grip, ASTM D2229 Adhesion (Steel tire cord & rubber) ****One Grip****

